

## Silent Feature

- Double and triple wall with stainless steel jacket 304 grade.
- Low water cut-off system
- User-friendly all-in-one type control panel
- Intelligent control system with PID program
- LCD display and soft touch operation key
- Over temperature Protection and auto-power break-off
- Digital timer

## Autoclave (Horizontal)



Model	SNS-HAC-100
Make	SN Solutions
Temperature Range	121°C to 135°C
Temperature Resolution	0.1°C
Wall Type	Double and Triple Wall
Working Temperature	121°C/135°C
Capacity	100 Liters
Dimension (Dia x Depth)	457.2 x 609.6 (mm)
internal & External Material	Stainless Steel 304 Grade
Controller	Microprocessor PID Controller
Sensor	PT-100
Display	LCD/LED
Design Type	Cylindrical