

Silent Feature

- Double and triple wall with stainless steel jacket 304 grade.
- Low water cut-off system
- User-friendly all-in-one type control panel
- Intelligent control system with PID program
- LCD display and soft touch operation key
- Over temperature Protection and auto-power break-off
- Digital timer

Autoclave (Horizontal)



Model	SNS-HAC-1200
Make	SN Solutions
Temperature Range	121°C to 135°C
Temperature Resolution	0.1°C
Wall Type	Double and Triple Wall
Working Temperature	121°C/135°C
Capacity	1200 Litters
Dimension (WxHxD)	914.4 x 914.4 x 1828.8 (mm)
internal & External Material	Stainless Steel 304 Grade
Controller	Microprocessor PID Controller
Sensor	PT-100
Display	LCD/LED
Design Type	Cylindrical