

Silent Feature

- **Controller temperature and RPM LCD display with timer and soft touch keypad Intelligent PID controlled**
- **Inside Illumination facility LED light**
- **Out Door Puf Insulated, fitted with silicon gasket and SS side lock and key**
- **Orbital Shaking Speed 30 to 250 RPM with various speed touch keypad**
- **Use only BOD**
- **Use only Shaker**

BOD Incubator Shaker



Model	SNS-BIS-175
Volume/Capacity	175
Removable Self	1
Temperature Range	5°C to 60°C
Temperature Accuracy	±0.2°C
Temperature Uniformity	±0.1°C at 80°C
External material	Mild steel powder coated
Internal Material	Stainless Steel 304 Grade
Motor	Sound less BLDC motor
Shaking	Sound less smooth shaking
Inner Door	Made up Acrylic Transparent glass
Shaking Capacity	16 x 250ml to 500 ml
Display	LCD
Power Supply	220-240V/50-60Hz Single Phase