

Silent Feature

- Made of Double walled Chamber-Inner made of Stainless Steel SS 304 grade and Outer made of Mild Steel powder coated with insulation of mineral wool/rock wool
- 1 shelves
- LCD Display with microprocessor PID controller with timer
- High and low temperature alarm
- Over temperature protection
- Air circulation fan

Digital Incubator



Model	SNS-BI-28
Make	SN Solutions
Capacity	28 ltrs
Inner Dimension	12 x 12 x 12 inch
Material of External Body	Mild Steel with Power Coated
Material of Internal Chamber	Stainless Steel 304 Grade
Material of Inner Door	Stainless Steel 304 Grade/Glass
Temperature range in °C	Ambient 60 °C.
Temperature Accuracy in °C	+/- 0.5°C
Temperature Uniformity	+/- 1 °C
Least Count	0.5 °C
Material of Insulation	Mineral Wool/Rock Wool