


Silent Feature

- Digital LCD display with back light .
- Big Size Ceramic hot plate with flame protection and durability .
- Heating and stirring can proceed simultaneously.
- Temperature & stirring speed is stepless adjusted.
- Real temperature & speed display on the screen simultaneously .
- Two sensors ,outside PT 100 sensor controls the liquid Temperature,inside K type sensor control the hot plate temperature (selectable)
- Over Temperature alarm
- **STANDARD ACCESSORIES**
 Stir Bar, Temperature Probe, Stand, Instruction Manual

Ceramic Hotplate Magnetic stirrer



Model	SNS-MS5LC
Make	SN Solutions
Display	Digital
Max. Volume	5000ml
Speed	100 - 2000rpm
Plate Size	190 x 190 mm
Top Plate Type	Ceramic
Max. Temp. (Plate Surface)	300°C
Temp. Accuracy	±1°C (<100°C) ±1°C (<300°C)
Power Wattage	600 W
Power supply	AC 220V / 50 Hz
Dimensions LxWxH (mm)	300 x 190 x 126
Weight	5 kg