

Silent Feature

- CO2 recovery rate At least of 6 min after door opening and closing event to attain 5% CO2
- Digital display for temperature CO2 and O2
- CO2 gas range 1 to 20%
- O2 gas range 1 to 19% with high quality O2 Sensor
- Incremental variation in CO2 gas range
- HEPA/ULPA filter for incoming gases
- Dry Heat Sterilization temperature in 180 degree Celsius
- Provision of display for temperature and alarm control
- Incubator should have auto-sterilization facility (Hot Air/Moist Air)
- CO2 cylinder and regulator provided with incubator

CO₂ Incubator



Model	SNS-CO-180
Volume/Capacity	180 Liters
Temperature Range	Ambient (+5to+60)
Display	LED Display
Temperature Control	Digital PID
Humidity Method	Natural vaporization >95%
Inner Chamber Material	Stainless Steel
Weight (Kg)	80
No. of Shelves	3
Co ² Range	0%~20%
Chamber Dimension (mm)(WxDxH)	473x528x710
Overall Dimensions (MM)	560x665x945
Power Supply	220-240V/50Hz Single Phase