

### Silent Feature

- 4 Inch attractive LCD Display
- CFC-free, HCFC-free refrigerents
- Auto tuning of controller
- Settable Hi/Low alarm points
- Power failure alarm
- Door open alarm
- Intelligent controller helps maintain temperature in case of sensor failure
- Histogram format of 24 hours temperature recording

### Deep Freezer



Model	SNS-DF-175
Volumu/Capacity	175 Liters
No. of Trays	4
No. of Internal Door	2
Cabinet type	Vertical
Dual compressor	YES
Operating Temperature Range	-40°C to -80°C
Temperature Control	Microprocessor with PT-100 Sensor
Display	4" LCD
Internal Body Material	SS 304(AISI Grade, Non Magnetic,)
External Body Material	Stainless Steel 304 Grade / Mild Steel Powder Coated
Noise Level	Less than 65 dB
Power Supply	220-240V/50Hz Single Phase