

Silent Feature

- Inner tank & outer body made of stainless steel
- Digital RED LED display
- Temperature with timer function
- Ceramic heaters provides improved heating effect
- Overheat protector, Cooling fan
- Moisture-proof PCB
- User friendly
- High-performance transducers

Digital Ultrasonic Cleaner



Model	SNS-UC-3
Make	SN Solutions
Tank Capacity	2.5 ltrs
Tank size (mm)	240 x 140 x 100
Overall size (mm)	270 x 170 x 240
Ultrasonic power	100 W
Heating power	100 W
Time setting	1-99 min.
Heating temperature	Ambient temperature to 80°C
Display	Digital RED LED
Frequency	40 kHz
Lid & Inner tank material	Stainless Steel 304 grade
Outer tank material	Stainless Steel 201 grade
Power supply	AC 220V/50Hz