

Silent Feature

- Double wall construction
- Double walled metal door with sponge - type silicon gasket for air-tight sealing
- Cord wire duly tested and inspected with stress factor as per CE standard
- Interior fabricated from high grade stainless steel (S.S.304).
- Exterior body made of Mild steel material - which is duly powder coated in attractive shades
- Full size polycarbonate sheet inner door to inspect samples without affecting the chamber temperature.

Environmental Chamber



| Model | SNS-EC-175 |
|-----------------------------------|--|
| Volumu/Capacity | 175 Liters |
| No. of Trays | 2 |
| Temperature Range | 5°C to 60°C |
| Humidity Range | 55% RH to 95% RH |
| Temperature measurement | PT 100 temperature sensor |
| Inner Chamber Size (W x D x H) | 50.5 x 41.5 x 83 (cm) |
| Display | LED / LCD / HMI Touch Screen |
| Inner Material | Stainless Steel 304 Grade |
| Outer Material | Mild Steel Powder Coated / Stainless Steel 304 Grade |
| Use only BOD | Yes |
| Controller | PID temperature controlled |
| No. of RH and Temperature Sensors | 1 |
| Power Supply | 220-240V/50Hz Single Phase |