

Silent Feature

- **Double wall construction**
- **Double walled metal door with sponge - type silicon gasket for air-tight sealing**
- **Cord wire duly tested and inspected with stress factor as per CE standard**
- **Interior fabricated from high grade stainless steel (S.S.304).**
- **Exterior body made of Mild steel material - which is duly powder coated in attractive shades**
- **Full size polycarbonate sheet inner door to inspect samples without affecting the chamber temperature.**

Environmental Chamber



Model	SNS-EC-280
Volumu/Capacity	280 Liters
No. of Trays	3
Temperature Range	5°C to 60°C
Humidity Range	55% RH to 95% RH
Temperature measurement	PT 100 temperature sensor
Inner Chamber Size (W x D x H)	57 x 55 x 87.5 (cm)
Display	LED / LCD / HMI Touch Screen
Inner Material	Stainless Steel 304 Grade
Outer Material	Mild Steel Powder Coated / Stainless Steel 304 Grade
Use only BOD	Yes
Controller	PID temperature controlled
No. of RH and Temperature Sensors	1
Power Supply	220-240V/50Hz Single Phase