

Silent Feature

- **Aerodynamic, Auto Bypass high constant volume type design**
- **Interior made of Chemical & Heat resistant, Fire Retardant, smooth finish, Easily Cleanable Panels made out of durable FRP/Phenolic resin construction integral work walls (6mm thick).**
- **Chemical resistant splash & spillage proof dished 'Jet Black Granite, worktop (18+ 1 mm thick) Skirting of 75 mm for all Side for no chemical spillage**
- **Vertical/Horizontal rising counter-balanced Toughened Float Glass sash (5 mm thick) smooth and light sash operation.**
- **20Ams MCB & 1 Set Relay & Contractor of reputed make**

Fume Hood



Model	SNS-FH-02
Make	SN Solutions
Design Type	Floor Model
Working Table	Granite Jet Black
Motor Power	1HP
Interior Illumination	1 Nos. 20 Watt LED Tube light Inside
Electrical Socket	4 Nos. with MCB (6/16A)
Inner Material	Stainless Steel 304 Grade
External Material	Stainless Steel 304 Grade / Mild Steel Powder Coteed / Wooden
Storage Cabinet	Yes
Openable Height	750 mm with Height Adjustment
Internal Dimension	790 x 600 x 870 mm approx
Outer Dimension	1000 x 800 x 2515 mm approx.
Power Supply	220V/50Hz