

**Silent Feature**

- **Aerodynamic, Auto Bypass high constant volume type design**
- **Interior made of Chemical & Heat resistant, Fire Retardant, smooth finish, Easily Cleanable Panels made out of durable FRP/Phenolic resin construction integral work walls (6mm thick).**
- **Chemical resistant splash & spillage proof dished 'Jet Black Granite, worktop (18+ 1 mm thick) Skirting of 75 mm for all Side for no chemical spillage**
- **Vertical/Horizontal rising counter-balanced Toughened Float Glass sash (5 mm thick) smooth and light sash operation.**
- **20Ams MCB & 1 Set Relay & Contractor of reputed make**

**Fume Hood**



<b>Model</b>	<b>SNS-FH-03</b>
<b>Make</b>	<b>SN Solutions</b>
<b>Design Type</b>	<b>Floor Model</b>
<b>Working Table</b>	<b>Granite Jet Black</b>
<b>Motor Power</b>	<b>1HP</b>
<b>Interior Illumination</b>	<b>1 Nos. 20 Watt LED Tube light Inside</b>
<b>Electrical Socket</b>	<b>4 Nos. with MCB (6/16A)</b>
<b>Inner Material</b>	<b>Stainless Steel 304 Grade</b>
<b>External Material</b>	<b>Stainless Steel 304 Grade / Mild Steel Powder Coteed / Wooden</b>
<b>Storage Cabinet</b>	<b>Yes</b>
<b>Openable Height</b>	<b>750 mm with Height Adjustment</b>
<b>Internal Dimension</b>	<b>990 x 600 x 870 mm approx</b>
<b>Outer Dimension</b>	<b>1200 x 800 x 2515 mm approx.</b>
<b>Power Supply</b>	<b>220V/50Hz</b>