

Silent Feature

- Caster Wheels
- Visibility of Item Stored
- Forst free Method of Defrosting
- Copper Cooling Coil Material
- Temprature Controler & ON/OFF Switch
- AC Rotating Compressor
- Refrigerants R404A
- 80mm Insulation Thickness
- 75 Minutes Holdover Time
- High Density PUF Insulation

Glass Door Refrigerator



Model	SNS-GDR-1100
Volumu/Capacity	1100 Litters
Shelves	10
Inner Door	Double
Position of Door	Front Opening
Number of Compressor	1
Operating Temperature Range	2°C ~ 8°C
Door Material	Glass Door
Temperature Control	Microprocessor PID with PT-100 Sensor
Type of cooling	Frost Free
Internal Body Material	Pre-coated mild steel
External Body Material	Prepainted Sheet
Dimension (WxDxH)	47 x 28 x 76
Power Supply	220-240V/50Hz Single Phase