

Silent Feature

- **Caster Wheels**
- **Visibility of Item Stored**
- **Forst free Method of Defrosting**
- **Copper Cooling Coil Material**
- **Temprature Controler & ON/OFF Switch**
- **AC Rotating Compressor**
- **Refrigerants R404A**
- **80mm Insulation Thickness**
- **75 Minutes Holdover Time**
- **High Density PUF Insulation**

Glass Door Refrigerator



Model	SNS-GDR-700
Volumu/Capacity	700 Litters
Shelves	08
Inner Door	Double
Position of Door	Front Opening
Number of Compressor	1
Operating Temperature Range	2°C ~ 8°C
Door Material	Glass Door
Temperature Control	Microprocessor PID with PT-100 Sensor
Type of cooling	Frost Free
Internal Body Material	Pre-coated mild steel
External Body Material	Prepainted Sheet
Dimension (WxDxH)	41 x 25 x 80
Power Supply	220-240V/50Hz Single Phase