

Silent Feature

- **Caster Wheels**
- **Visibility of Item Stored**
- **Forst free Method of Defrosting**
- **Copper Cooling Coil Material**
- **Temprature Controler & ON/OFF Switch**
- **AC Rotating Compressor**
- **Refrigerants R404A**
- **80mm Insulation Thickness**
- **75 Minutes Holdover Time**

Glass Door Refrigerator



| Model | SNS-GDR-200 |
|-----------------------------|----------------------------|
| Volumu/Capacity | 200 Litters |
| Shelves | 3 |
| Inner Door | Single |
| Position of Door | Front Opening |
| Number of Compressor | 1 |
| Operating Temperature Range | 2°C~8°C |
| Door Material | Glass Door |
| Temperature Display | yes |
| Microprocessor Controller | yes |
| Internal Body Material | Pre-coated mild steel |
| Lock | Yes |
| Dimension (WXDXH) | 22*22*65 |
| Power Supply | 220-240V/50Hz Single Phase |