

Silent Feature

- **Caster Wheels**
- **Visibility of Item Stored**
- **Forst free Method of Defrosting**
- **Copper Cooling Coil Material**
- **Temprature Controler & ON/OFF Switch**
- **AC Rotating Compressor**
- **Refrigerants R404A**
- **80mm Insulation Thickness**
- **75 Minutes Holdover Time**
- **High Density PUF Insulation**

Glass Door Refrigerator



Model	SNS-GDR-400
Volumu/Capacity	400 Litters
Shelves	5
Inner Door	Single
Position of Door	Front Opening
Number of Compressor	1
Operating Temperature Range	2°C~8°C
Door Material	Glass Door
Temperature Display	yes
Microprocessor Controller	yes
Internal Body Material	Pre-coated mild steel
Lock	Yes
Dimension (WXDXH)	25 x 24 x 76
Power Supply	220-240V/50Hz Single Phase