

Silent Feature

- Digital display of temperature
- Vacuum display in for 3 colors red, yellow and green. Red color show is not vacuum, yellow color show is not proper vacuum and green color show is fine vacuum
- CFC and HCFC Refreigerant
- Exhaust oil mist filter
- 8-port stainless steel manifold
- High water vapor handling capacity
- Condenser material Stainless steel 316 grade

Lyophilizer (Freeze Dryer)



Model	SNS-DF-40
Material Capacity	3.5 ltrs.
Ice Removal Capacity in 24 Hrs	2.3 liters
Condenser Temperature	-55°C
Vacuum & Temp. Controller	Microprocessor Based Digital Controller
Material of Chamber	SS 316 Grade
Outer material	Mild steel powder coated
Foot space	650 x 850mm
Max. Ice capacity	1.5 kg
Refrigerant	CFC free and PUF insulated
Vacuum & Temp. Indications	Digital
Vacuum Pump Capacity	250 L/Min
Power Supply	220-230V/50Hz Single Phase