


Silent Feature

- Manual Temperature Compensation.
- There is a microprocessor chip inside the meter and the appearance of meter looks nice with user friendly operation. Some of the features are as under.
- Large LCD display with blue background simultaneous display of pH AND temperature.
- It has intelligent functions for calibration, data storage / Memory.
- Automatically display the percentage of slope after calibration is over .
- Equipped with electrode holder and calibration buffer solution, which is more convenient to use.

pH Meter



Model	SNS-pH-09
Make	SN Solutions
pH Range	0 to 14.00 pH
pH Resolution	0.1 ph / 0.01 pH
Accuracy	±0.01 pH
mV Range	±1999 mV
Calibration Points	up to 3 points
mV Resolution	1 mV
Temperature Range	0 to 100°C
Temperature Accuracy/Resolution	±0.5°C/0.1°C
Display	LCD (135 x 75 mm)
Power Requirements	DC9V, using AC adapters, AC 220V / 50Hz
Dimension (L x W x H)	160 x 190 x 170 (mm)