

Silent Feature

- Large backlight LCD display show mode indicator and controls are required.
- Up to 3 points push-button calibration with auto buffer recognition.
- Selectable pH buffer standards and temperature units.
- Automatically calculates and displays electrode slope after calibration.
- Automatic temperature compensation ensures highly accurate measurements over the entire range.
- Automatic endpoint function freezes the stale measure value.
- Setup menu lets user customizes buffer standards, calibration, temperature units, automatic and manual hold function, auto-off feature to meet personal preferences.
- Reset function can resume all setting to factory default values

pH Meter



Model	SNS-pH-10
Make	SN Solutions
pH Range	-2.00 to 14.00 pH
pH Resolution	0.01 pH
Accuracy	±0.01 pH
mV Range	±1999 mV
mV Accuracy	±1 mV
mV Resolution	1 mV
Temperature Range	0 to 105°C
Temperature Accuracy/Resolution	±0.5°C/0.1°C
Display	LCD (135 x 75 mm)
Hold & Auto Function	Available
Dimension (L x W x H)	210 x 205 x 75 (mm)