

Silent Feature

- Microprocessor PID Control / Auto tubing / Calibration.
- User-friendly door system
- Air flow system very efficient for plant growth
- Heat rise restraint
- Easy cleaning of condenser air filter
- Over current protection.
- Door opening alarm
- Low water level alarm
- Air circulation fan for smooth airflow

Plant Growth Chamber



Model	SNS-PGC-450
Volume/Capacity	450 Litters
Digital timer	1 min. to 99 hr. 59 min
Temperature	Fluctuation: $\pm 0.5^{\circ}\text{C}$, Variation: $\pm 0.8^{\circ}\text{C}$
Illumination Range	0 to 20,000 Lux
Control	24 Program
Lamp	Led light / fluorescent Light
Controller Type	LCD/LED/ HMI/PLC controller with graph
Display	HMI Touch Screen Display
Humidity Range	60-95% RH at 30-50- $^{\circ}\text{C}$
Compressor	1/2HP
Temperature Heater power	500Wx2EA
Humidity Heater power	1300W