

Silent Feature

- PT 100 temperature sensor
- Both side folding handle for easy movement
- overheat protection and liquid level alarm function
- stainless steel inner tank
- Water level function
- Fast cooling, Low noise, Low vibration, Stable & Reliable performance
- LED, Microprocessor, high precision, PID temperature control system

Refrigerated Water Bath



Model	SNS-RWB-12
Make	SN Solutions
Capacity	12.0 ltrs
Temperature Range	5°C to 90°C
Material of Internal Body	Stainless Steel 304 Grade
Material of External Body	Mild steel duly epoxy powder coated
Internal Dimension (H x W x D)	240 x 250 x 200 mm
Accuracy	± 0.5°C to 50°C
Display	LCD Display
Pump Circulation	Internal & External
Power Supply	220 to 240V AC
Load	2.0 kw