

Silent Feature

- Inner Material Stainless Steel 304 Grade
- Outer Material Mild Steel Powder Coated
- Single or Double Solid Doors with Lock & Key
- Acrylic Door
- CFC Free Refrigerant
- Digital PID Controller for Temperature & Humidity System
- 24x7 weekly digital timer
- Kanthal A-1 Tubular Heater
- Audio/Visual Alarm for Temperature and Humidity Deviation

Seed Germinator



Model	SNS-SG-340
Make	SN Solutions
Capacity	340 Itrs.
Temperature Range	5°C to 60°C
Temperature and Humidity Controller	Digital PID controller
Temperature Fluctuation	±1°C
Temperature Sensor	PT-100
Humidity Range	50% to 95% RH
Insulation	75 mm PUF
Display	LED Display
Light	Fluorescent Lamp
Shelves	4
Dimension (WxDxH)	605 x 555 x 910 mm
Power	220-240 V, 50/60 Hz, Single Phase

