

Silent Feature

- User Password Protection available.
- Microprocessor controlled and completely programmable.
- Large LCD display with Time setting.
- Temperature Indicator and controller.
- Auto-tuning for convenience and optimal processing efficiency.
- Variable Power Output control.
- Auto frequency tracking.
- Output amplitude 0-100% adjustment.
- Automatic fault alarm.

Ultrasonic Probe Sonicator



Model	SNS-PRO-1200
Make	SN Solutions
Frequency	19.5-20.5KHz
Power	12 – 1200W Adjustable
Crushing Capacity	3-1000 ml
Standard Probe	20 mm
Optional Probes	10/15/22 mm
Display	7" Colorful Touch Screen Display
Temperature Range	0-100°C
Waveform Display	Pulse, Frequency and Time
Data Storage	20 sets
Dimension (L*W*H)	470*305*430 mm

