

Silent Feature

- User Password Protection available.
- Microprocessor controlled and completely programmable.
- Large LCD display with Time setting.
- Temperature Indicator and controller.
- Auto-tuning for convenience and optimal processing efficiency.
- Variable Power Output control.
- Auto frequency tracking.
- Output amplitude 0-100% adjustment.
- Automatic fault alarm.

Ultrasonic Probe Sonicator



| Model | SNS-PRO-650 |
|-------------------|----------------------------------|
| Make | SN Solutions |
| Frequency | 20-25KHz |
| Power | 6.5 – 650W Adjustable |
| Crushing Capacity | 0.2-500 ml |
| Standard Probe | 6 mm |
| Optional Probes | 2/3/6/8/10/12 mm |
| Display | 7" Colorful Touch Screen Display |
| Temperature Range | 0-100°C |
| Waveform Display | Pulse, Frequency and Time |
| Data Storage | 20 sets |
| Dimension (L*W*H) | 470*305*430 mm |

