

Vacuum Oven



Model	SNS-VO-45
Make	SN Solutions
Capacity	45 ltrs.
Finish for Outer Housing Material	Stainless Steel 304 Grade
Outer Housing Material	Mild Steel
Vacuum Chamber Material	Stainless Steel
Dimension (L x H x W)	355 x 355 x 355 mm
Maximum Operating Temperature	50°C
Minimum Operating Temperature	-10°C
Shelves	1
Controller Type	Digital PID Controller
Ingress Protection	Ip54
Power Supply	220V/50Hz