

Silent Feature

- Double Wall with Stainless Steel Jacket 304 Grade.
- Triple Safety Feature
- HMI Display with Touch Screen
- Low Water Cut-off System
- User Friendly all-in-one type Control Panel
- Intelligent Control System with PID Program
- Open the Door Heater off Automatically
- Over Temperature Protection and Auto-power Break-off
- Chamber Temperature, Alarm Indicator, Low Water Indicator
- Inbuilt MCB in the Device
- Temperature Graph, Data Logger and USB Port

Autoclave (Vertical)



Model	SNS-VAC-22
Make	SN Solutions
Capacity	22 Liters
Temperature Range	121°C to 135°C
Temperature Resolution	0.1°C
Temperature Uniformity	1°C
Working Temperature	125°C
Range of Pressure Gauge	30.0 (Pound per square inch)
Dimension (Dia x Depth)	250x450 (mm)
Inner Chamber Material	Stainless Steel 304 Grade
Outer Chamber Material	Stainless Steel 304 Grade
Digital Timer	1-99 minutes
Sensor	PT-100
Electric Load	1.5 KW