

## Silent Feature

- Double Wall with Stainless Steel Jacket 304 Grade.
- Triple Safety Feature
- HMI Display with Touch Screen
- Low Water Cut-off System
- User Friendly all-in-one type Control Panel
- Intelligent Control System with PID Program
- Open the Door Heater off Automatically
- Over Temperature Protection and Auto-power Break-off
- Chamber Temperature, Alarm Indicator, Low Water Indicator
- Inbuilt MCB in the Device
- Temperature Graph, Data Logger and USB Port

## Autoclave (Vertical)



| Model                   | SNS-VAC-78                   |
|-------------------------|------------------------------|
| Make                    | SN Solutions                 |
| Capacity                | 78 Liters                    |
| Temperature Range       | 121°C to 135°C               |
| Temperature Resolution  | 0.1°C                        |
| Temperature Uniformity  | 1°C                          |
| Working Temperature     | 125°C                        |
| Range of Pressure Gauge | 30.0 (Pound per square inch) |
| Dimension (Dia x Depth) | 400x600 (mm)                 |
| Inner Chamber Material  | Stainless Steel 304 Grade    |
| Outer Chamber Material  | Stainless Steel 304 Grade    |
| Digital Timer           | 1-99 minutes                 |
| Sensor                  | PT-100                       |
| Electric Load           | 3 KW                         |