

Silent Feature

- **MULTIPLE RESULTS READOUT** - Can display wavelength, absorption and transmittance with 5 results per screen. It also has a memory store of up to 200 results.
- **AUTO SETTING WAVELENGTH** - Users set wavelength automatically through arrow keys to avoid operation errors.
- **STANDARD CURVE** - Can set up various standard curves according to customer's solutions and find the concentration of unknown solutions
- **DATA OUTPUT** - SNS-SPV325 is equipped with USB port to connect with a PC to display spectrum scanning, kinetics and multi wavelength testing results on the screen.
- **2.5" LCD SCREEN** - Equipped with a 2.5 LCD screen to give a clear display of standard curves and groups of results.

Visible Spectrophotometer



Model	SNS-SPV325
Make	SN Solutions
Cell Holder	4-position 10mm cell holder
Bandwidth	4nm
Optical System	Single beam, Grating 1200 lines/mm
Wavelength Range	325-1000nm
Wavelength Accuracy	±1nm
Wavelength Repeatability	±0.5nm
Dimension (LxWxH)	420x280x180 mm
Photometric Accuracy/ Repeatability	±0.5%T / ±0.3%T
Data Output Port	USB Port & Parallel Port(Printer)
Lamps	Tungsten Lamp
Detector	Silicon Photodiode
Power supply	AC 220V / 50Hz