

Silent Feature

- A double walled bath with inner chamber made of SS sheet (304 grade) Argon Arc Welded, and outer wall made of MS sheet duly epoxy powder coated.
- Provided with SS lid having 75mm dia.
- Holes with removable concentric rings.
- Provided with on/off switch, indicator, power cord etc.
- The gap between outer and inner walls is filled with special grade glass wool insulation.
- Temperature from +5°C ambient to 99°C.

Water Bath Rectangular (Double Wall)



Model	SNS-WBR-06
Make	SN Solutions
Capacity	08 ltrs.
Temperature Range	+5°C Ambient to 99°C
Internal Material	Stainless Steel 304 Grade
External Material	Stainless Steel 304 Grade / Mild Steel Powder Coteed
No. of Holes	6 Nos.
Hole Diameter	75 mm
Heating Element	1.0 KW
Controller	Microprocessor PID
Display	LED Display
Power supply	220V. 50 Hz single phase AC