

Silent Feature

- A double walled bath with inner chamber made of SS sheet (304 grade) Argon Arc Welded, and outer wall made of MS sheet duly epoxy powder coated.
- Provided with SS lid having 75mm dia.
- Holes with removable concentric rings.
- Provided with on/off switch, indicator, power cord etc.
- The gap between outer and inner walls is filled with special grade glass wool insulation.
- Temperature from +5°C ambient to 99°C.

Water Bath Rectangular (Double Wall)



Model	SNS-WBR-12
Make	SN Solutions
Capacity	15 ltrs.
Temperature Range	+5°C Ambient to 99°C
Internal Material	Stainless Steel 304 Grade
External Material	Stainless Steel 304 Grade / Mild Steel Powder Coteed
No. of Holes	12 Nos.
Hole Diameter	75 mm
Heating Element	2.0 KW
Controller	Microprocessor PID
Display	LED Display
Power supply	220V. 50 Hz single phase AC