

Silent Feature

- Properly insulated. Heating elements are tubular type, made of copper duly nickel plated.
- Covered with perforated SS Tray.
- Temperature is controlled with Temperature controller is graduated in degree centigrade.
- Provided with on/off switch, indicator, power cord etc.
- Temperature Controller Led Display

Water Bath Shaker



Model	SNS-WBS-09
Make	SN Solutions
Temperature Range	Ambient to 100°C
Internal Material	Stainless Steel 304 Grade
External Material	Stainless Steel 304 Grade / Mild Steel Powder Coteed
Flask	9
Shaking Capacity	250/500ml
Drain Valve	Yes (Inbuilt)
Controller	Microprocessor PID
Display	LED/LCD Display
Power supply	220V. 50 Hz single phase AC