

Silent Feature

- Properly insulated. Heating elements are tubular type, made of copper duly nickel plated.
- Covered with perforated SS Tray.
- Temperature is controlled with Temperature controller is graduated in degree centigrade.
- Provided with on/off switch, indicator, power cord etc.
- Temperature Controller Led Display
- Double walled with outer body made up of PCGI sheet duly powder coated and inner chamber and lid is made up of stainless steel sheet (SS 304 grade) .

Water Bath Shaker



| Model | SNS-WBS-12 |
|-------------------|--|
| Make | SN Solutions |
| Temperature Range | Ambient to 100°C |
| Internal Material | Stainless Steel 304 Grade |
| External Material | Stainless Steel 304 Grade / Mild Steel Powder Coteed |
| Flask | 12 |
| Shaking Capacity | 250/500ml |
| Drain Valve | Yes (Inbuilt) |
| Controller | Microprocessor PID |
| Display | LED/LCD Display |
| Power supply | 220V. 50 Hz single phase AC |