

Silent Feature

- PT100 temperature sensor
- Both sides folding handle for easy movement
- Overheat protection and liquid level alarm function.
- Stainless steel inner tank, Water level alarm function
- Fast cooling, low noise, low vibration, stable and reliable performance
- LCD, Microprocessor, High Precision, PID temperature control system.

Recirculating Chiller



Model	SNS-RC-06
Make	SN Solutions
Tank Capacity	06 Liters
Temperature Range	-20°C TO 100°C
Temperature Stability	±1°C
Pump circulation	Internal & External
Temperature sensor	Pt100
Display	LCD
Temperature Control	Microprocessor PID
Refrigeration	Air Cooled refrigeration through compressor
Tank Size (L x W x D)	180x 290x 150mm
Bath dept	150 mm
Weight	44Kg