

## Silent Feature

- PT100 temperature sensor
- Both sides folding handle for easy movement
- Stainless steel inner tank, Waterlevel alarm function.
- Fast cooling, low noise, low vibration, stable and reliable performance.
- LCD, Microprocessor, High Precision, PID temperature control system

## Recirculating Chiller



Model	SNS-RC-12
Make	SN Solutions
Tank Capacity	12 Liters
Temperature Range	-20°C TO 100°C
Temperature Stability	±1°C
Pump circulation	Internal & External
Temperature sensor	Pt100
Display	LCD
Temperature Control	Microprocessor PID
Refrigeration	Air Cooled refrigeration through compressor
Tank Size (L x W x D)	240 x 250x 200mm
Bath dept	200mm
Weight	50 Kg