

Silent Feature

- **PT100 temperature sensor**
- **Both sides folding handle for easy movement**
- **Stainless steel inner tank, Waterlevel alarm function.**
- **Fast cooling, low noise, low vibration, stable and reliable performance.**
- **LCD, Microprocessor, High Precision, PID temperature control system**

Recirculating Chiller



Model	SNS-RC-14
Make	SN Solutions
Tank Capacity	14 Liters
Temperature Range	1°C TO 100°C
Temperature Stability	±1°C
Pump circulation	Internal & External
Temperature sensor	Pt100
Display	LED Display
Temperature Control	Microprocessor PID
Refrigeration	Air Cooled refrigeration through compressor
Tank Size (L x W x D)	270 x 240x 240mm
Bath dept	200mm
Weight	60 Kg