

## Silent Feature

- Intelligent Microprocessor PID Controlled and PT-100 Sensor
- LED/LCD/HMI Touch Screen Display (Selectable)
- Automatic Visual & Audible alarm for Over Temperature
- Out door Polyurethane/PUF Insulated, Fitted with Silicon Gasket and Stainless Steel Side Lock and Key
- Inside Illumination Fluorescent Lamp / LED light facility
- Manual/Automatic Defrost System

## BOD Incubator



Model	SNS-BOD-112	SNS-BOD-175	SNS-BOD-220	SNS-BOD-280	SNS-BOD-340
Volume/Capacity	112 ltrs./4 (Cu.ft)	175 ltrs./6 (Cu.ft)	220 ltrs./8 (Cu.ft)	280 ltrs./10 (Cu.ft)	340 ltrs./12 (Cu.ft)
Temperature Range	+5°C to 60°C				
Temp. Uniformity	± 0.5°C				
Temp. Accuracy	± 0.5°C				
Volume/Capacity	2 Trays	2 Trays	3 Trays	3 Trays	3 Trays
Air Circulation	Motorized Blower				
Exterior Body Material	Mild Steel Powder Coated / Stainless Steel Powder Coated (Selectable)				
Exterior Body Material	Mild Steel Powder Coated / Stainless Steel Powder Coated (Selectable)				
Insulation Material	Polyurethane / PUF (Selectable)				
Outer Door Material	Mild Steel Powder Coated / Stainless Steel Powder Coated (Selectable)				
Display	LED / LCD / HMI Touch Screen (Selectable)				
Power Supply	220 - 240 V / 50 Hz Single Phase				