

Silent Feature

- High-quality ISI Mark Stainless steel 304 Grade .
- Gas and vacuum cock
- Convenient for cleaning.
- Fluorescent tubes Illumination 1000 Lumens.
- Air pressure gauge.
- Noise Level 60 decibels.
- LCD display
- Short circuit breaker .
- Provided with two pre-filters.

Biosafety Cabinet



Model	SNS-BSC-04
Make	SN Solutions
Internal Material	Stainless Steel 304 Grade
Filtration of particulate	0.3 micron
Main Filter	HEPA Filter 99.999%Efficient
Exhaust Filter	HEPA Filter 99.999%Efficient
Duct Size	6 inch
Virus Protraction	Virus burn heater
Working size	4 x 2 x 2 (ft)
Outer size	4.6 x 2.9 x 6 (ft)
Fan Motor Type	Direct Drive Blower Motor
Power Supply	220-230V/50-60Hz
Gross Weight	110 Kg