

### Silent Feature

- **Made of Double Walled Chamber-Inner made of Stainless Steel SS 304 Grade and Outer made of Mild Steel Powder Coated with Insulation of Glass Wool**
- **1 Shelves**
- **LCD Display with Microprocessor PID Controller with Timer**
- **High and Low Temperature Alarm**
- **Over Temperature Protection**
- **Air Circulation Fan**
- **Show in Display Set Value (SV) and Present Value (PV)**
- **Time , Temperature Setting and Temperature Reading**

### Hot Air Oven



Model	SNS-HO- 28
Make	SN Solutions
Capacity	28 ltrs
Inner Dimension	12 x 12 x 12 inch
Material of External Body	Mild Steel with Power Coated
Material of Internal Chamber	Stainless Steel 304 Grade
Material of Inner Door	Stainless Steel 304 Grade/Glass
Temperature range in °C	Ambient to +5°C to 250 °C.
Temperature Accuracy in °C	±1°C
Air Circulation	Motorized Blower for Force Air Circulation
Temperature Sensor	PT 100
Power supply	220V, 50Hz Single Phase
Material of Insulation	Glass Wool