

Silent Feature

- Made of Double Walled Chamber-Inner made of Stainless Steel SS 304 Grade and Outer made of Mild Steel Powder Coated with Insulation of Glass Wool
- 2 Shelves
- LCD Display with Microprocessor PID Controller with Timer
- High and Low Temperature Alarm
- Over Temperature Protection
- Air Circulation Fan
- Show in Display Set Value (SV) and Present Value (PV)
- Time , Temperature Setting and Temperature Reading

Hot Air Oven



Model	SNS-HO-95
Make	SN Solutions
Capacity	95 ltrs
Inner Dimension	18x 18 x 18 inch
Material of External Body	Mild Steel with Power Coated
Material of Internal Chamber	Stainless Steel 304 Grade
Material of Inner Door	Stainless Steel 304 Grade/Glass
Temperature range in °C	Ambient to +5°C to 250 °C.
Temperature Accuracy in °C	±1°C
Air Circulation	Motorized Blower for Force Air Circulation
Temperature Sensor	PT 100
Power supply	220V, 50Hz Single Phase
Material of Insulation	Glass Wool