

Silent Feature

- LCD Display
- Soundless Motor
- RPM Push Key
- Start Stop Key
- Timer
- On/Off Switch
- Autorun

Orbital Shaker



Model	SNS-RFS-16
Make	SN Solutions
Outer Material	Heavy Duty MS Angle Iron Frame
Motor Power	0.5 HP
RPM	30 to 250
RPM Accuracy	±2 RPM
Timer	1 to 999 minutes
Motor Type	Brushless DC Motor
Inner Size	17 x 17 Inch
Soft Key	RPM, Timer, Start and Stop key
Flask Capacity	25 x 250ml to 500 ml
Power Supply	220V/50Hz