

Silent Feature

- **Made of Double Walled Chamber-Inner made of Stainless Steel SS 304 Grade and Outer made of Mild Steel Powder Coated with Insulation of Glass Wool**
- **4 Shelves**
- **LCD Display with Microprocessor PID Controller with Timer**
- **High and Low Temperature Alarm**
- **Over Temperature Protection**
- **Air Circulation Fan**
- **Show in Display Set Value (SV) and Present Value (PV)**
- **Time , Temperature Setting and Temperature Reading**

Hot Air Oven



| Model | SNS-HO- 335 |
|------------------------------|--|
| Make | SN Solutions |
| Capacity | 335 ltrs |
| Inner Dimension | 24 x 24 x 36 inch |
| Material of External Body | Mild Steel with Power Coated |
| Material of Internal Chamber | Stainless Steel 304 Grade |
| Material of Inner Door | Stainless Steel 304 Grade/Glass |
| Temperature range in °C | Ambient to +5°C to 250 °C. |
| Temperature Accuracy in °C | ±1°C |
| Air Circulation | Motorized Blower for Force Air Circulation |
| Temperature Sensor | PT 100 |
| Power supply | 220V, 50Hz Single Phase |
| Material of Insulation | Glass Wool |