

### Silent Feature

- Magenagalic Gauge and UV Hour Meter
- Work Zone Stainless Steel Working Area
- Window and Side Wall Material Acrylic Transparent
- LCD Display
- LED Light Intensity 800 lux
- Mini plate HEPA
- DOP HEPA filter testing
- Stainless Steel 304 grade working area
- Digital Pressure Gauge with Alarm
- Vertical sliding type Front Door
- Noise level ( $\pm 5$ ) in 65dB
- Door and Frames hermatically sealed

### Vertical Laminar Air Flow



Model	SNS-LAF-04
Make	SN Solutions
Type of Configuration	Vertical
Fan Motor Type	Direct Drive Blower Motor
Cabinet Construction	Stainless Steel 304 Grade
Working table	Stainless Steel 304 Grade
Main Filter	HEPA Filter 99.995% Efficient
Pre Filter	Efficiency of 90% down to 0.3 Micron size
LED Light	30 W, 1 Pcs
Side wall	Acrylic Transparent
Working Size (L x W x H)	4 x 2 x 2 (ft)
Window material	Acrylic Transparent
Power Supply	100-240VAC/50-60Hz
Gross Weight	80 Kg